

www.flamor.com



NECKLITE

moldable neck brace

Introducing a moldable emergency neck brace, comprised of soft aluminum and coated with high density PU foam. These specific design characteristics allow the brace to be folded or rolled for compact storage. This allows for convenient transport in backpacks, kits and narrow pockets. A series of moldable small wings are positioned under the chin in order to conform to each patient's unique features. Once applied to a patient, it becomes rigid and the neck is tightly immobilized.

**ONE SIZE
NO SETTINGS
FIT & MOULD**



Cargo Pocket



Small bag kit



Backpacks



Vest



**COMPACT - LIGHT
FAST - EFFECTIVE
X-RAY TRANSLUCENT
EASY TO USE**



Thin



Moldable



- The device is light and thin, measuring 5 mm thick with a height of 180 mm and a total length of 550 mm.
- The locking system uses a 5 cm black Velcro strap.
- Once in place, the device features an elliptical hole on the backside and a large oval shaped hole on the front side which allows for a tracheotomy.



X-Ray translucent



Easy to locate



Military colors:

Outside: Dark grey (low visibility)

Inside: Dark grey (low visibility)

Civilian colors:

Outside: Yellow - **Inside:** Dark grey


Weight and dimensions:

Weight: 105 gr.

Width: 550 mm

Height: 180 mm

Thickness: 5 mm



Pediatric colors:

Outside: Blue

Inside: Dark grey


Weight and dimensions:

Weight: 75 gr.


Width: 440 mm

Height: 144 mm


Thickness: 5 mm



Slide the back of the collar, with the elliptical hole, behind the patient's neck. Continue pushing beneath the patient until you can reach it on the opposite side.



Wrap the front side around the lower jaw, aligning the oval hole under the chin. Secure the collar tightly around the patient's neck using the Velcro strap.



Push the small wings against the chin and check that they are in a supportive position. Adjust the collar into a position that immobilizes the neck of the patient and readjust the Velcro strap if necessary.



Via Biandrate, 80 - 28060 SAN PIETRO MOZZO (NO) - ITALY
Tel. +39 0321 442441 - Fax +39 0321 392190 - E-mail: info@flamor.com